




Email : sales@laboratoryinstrumentindia.com

Product Name :
Automotive Measurement Multimeter Trainer

Product Code :
LBNY-0005-120000194



Description :

Features:

Six individual errors can be selected by rotary switches in each circuit
For troubleshooting and voltage drop measurements 47 test sockets are available
Automotive relation with three typical, simple circuits for automotive lamps: conventional, with relay, with transistor
Colored printed circuits in standard-compliant display

Learning Objectives:

Use part voltage measurements in intact and defective circuits as a diagnostic method
Measuring in the correct place in the harness
Creating and evaluating protocols for reference values and actual values
Get to know typical lamp circuits
Using the multimeter to measure voltage
Interpret measured voltages correctly and thus identify the main cause of an error.

Technical Specification :

Dimensions: L x W x H 297 x 266 x 80 mm
Power supply: 12 Volt DC
Weight: about 2 kg.



Laboratory Instrument India