




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Automotive Measurement Multimeter Trainer

**Product Code :**  
LBNY-0005-120000194



#### **Description :**

##### **Features:**

Six individual errors can be selected by rotary switches in each circuit  
For troubleshooting and voltage drop measurements 47 test sockets are available  
Automotive relation with three typical, simple circuits for automotive lamps: conventional, with relay, with transistor  
Colored printed circuits in standard-compliant display

##### **Learning Objectives:**

Use part voltage measurements in intact and defective circuits as a diagnostic method  
Measuring in the correct place in the harness  
Creating and evaluating protocols for reference values and actual values  
Get to know typical lamp circuits  
Using the multimeter to measure voltage  
Interpret measured voltages correctly and thus identify the main cause of an error.

#### **Technical Specification :**

Dimensions: L x W x H 297 x 266 x 80 mm  
Power supply: 12 Volt DC  
Weight: about 2 kg.



**Laboratory Instrument India**