



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Automotive DC Converter Trainer

**Product Code :**  
LBNY-0005-120000198



**Description :**

LCDs to display the respective input and output voltages  
Adjustable load resistors in both circuits  
Voltage measurement points upstream and downstream of all components, including additional current-measuring points on all components  
DC/DC voltage converters are the link between high-voltage batteries and three-phase machines. The innovative laboratory trainer for step-down and step-up conversion enables the detailed study of this interaction.  
Fine adjustable signal generators  
Low voltage for safe handling  
Colored printed, clear presentation of the two circuits to show function and structure of DC/DC converters.

**Technical Specification :**

Dimensions: L x W x H 297 x 266 x 80 mm  
Power supply: 12 Volt DC  
Weight: about 1 kg.



**Laboratory Instrument India**

