



Email : sales@laboratoryinstrumentindia.com

Product Name :
Automotive DC Converter Trainer

Product Code :
LBNY-0005-120000198



Description :

LCDs to display the respective input and output voltages

Adjustable load resistors in both circuits

Voltage measurement points upstream and downstream of all components, including additional current-measuring points on all components

DC/DC voltage converters are the link between high-voltage batteries and three-phase machines.

The innovative laboratory trainer for step-down and step-up conversion enables the detailed study of this interaction.

Fine adjustable signal generators

Low voltage for safe handling

Colored printed, clear presentation of the two circuits to show function and structure of DC/DC converters.

Technical Specification :

Dimensions: L x W x H 297 x 266 x 80 mm

Power supply: 12 Volt DC

Weight: about 1 kg.



Laboratory Instrument India

