



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Car Electronics Trainer

**Product Code :**  
LBNY-0005-120000199



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

A prerequisite for secure diagnosis is knowledge of the component Specifications and competence to verify its function by measurements. Information, measuring circuits and assessments focus on vehicle-related situations out of the great variety of application options. The system comprises all required components it is compact, can be used at once, and permits free wiring. The focus of the pedagogically valuable course documents is on car diagnosis. In vehicle systems, very different electronics components are integrated in circuits. The solutions book contains additional further explanations and solution charts. The unique course comprises comprehensive documentation and extensive worksheets for teachers and students.

**Technical Specification :**

Dimensions:  
Width 300 mm  
Depth 265 mm  
Height approx. 80 mm  
Weight approx. 1 kg



**Laboratory Instrument India**

