



Email : sales@laboratoryinstrumentindia.com

Product Name :
Automotive Oscilloscope Trainer

Product Code :
LBNY-0005-120000200



Description :

Features:

Six individual errors can be selected by rotary switch for both components
Sine, square, and pulse width modulated control of components
Colored printed circuits in standard-compliant display
Automotive relating with two classic components: heating element and electric motor
For troubleshooting and measurement 35 test sockets are available

Learning Objectives:

Using the oscilloscope to measure voltage
Measuring in the correct place in the harness
With the help of sine and square wave signals you learn the fundamentals for effective value (RMS) and learn how to use your gauge
Insert part voltage measurements in intact and defective circuits as a diagnostic method
Rate measured signals correctly and thus identify the main cause of an error.

Technical Specification :

Power supply: 110/220 Volt AC
Dimensions: L x W x H 297 x 266 x 80 mm
Weight: about 3 kg



Laboratory Instrument India