



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Automotive Oscilloscope Trainer

**Product Code :**  
LBNY-0005-120000200



#### **Description :**

##### **Features:**

Six individual errors can be selected by rotary switch for both components  
Sine, square, and pulse width modulated control of components  
Colored printed circuits in standard-compliant display  
Automotive relating with two classic components: heating element and electric motor  
For troubleshooting and measurement 35 test sockets are available

##### **Learning Objectives:**

Using the oscilloscope to measure voltage  
Measuring in the correct place in the harness  
With the help of sine and square wave signals you learn the fundamentals for effective value (RMS) and learn how to use your gauge  
Insert part voltage measurements in intact and defective circuits as a diagnostic method  
Rate measured signals correctly and thus identify the main cause of an error.

#### **Technical Specification :**

Power supply: 110/220 Volt AC  
Dimensions: L x W x H 297 x 266 x 80 mm  
Weight: about 3 kg



**Laboratory Instrument India**