



Email : sales@laboratoryinstrumentindia.com

Product Name :
Power Steering Trainers

Product Code :
LBNY-0005-120000210



www.laboratoryinstrumentindia.com

Description :

An ampere meter provides an indication of the load and an in-line hydraulic pressure gauge monitors the hydraulic fluid circuit. Produces two types of Power Steering Trainers to demonstrate the principles and operation of power steering systems. In both trainers an electrical motor drives the hydraulic pump with belt drive. A resistance mechanism provides realistic operation and gives a load to the hydraulic circuit.

Technical Specification :

Power Steering Trainers



Laboratory Instrument India