



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Brake Chassis Trainer

Product Code :
LBNY-0005-120000211



Description :

The chassis includes an integrated floor stand that is mounted on wheels for mobility. Hydraulic Brake Chassis Trainer is designed for student skills training in the operation, installation, maintenance and repair of hydraulic brake systems. A color schematic chart of the brake system and circuit is mounted onto the trainer for reference. The unit consists of original automotive components which are mounted into a simulated model auto chassis.

Technical Specification :

Hydraulic Brake Chassis Trainer



Laboratory Instrument India