



Email : sales@laboratoryinstrumentindia.com

Product Name :
Shock Absorber Test Stand

Product Code :
LBNY-0005-120000212



www.laboratoryinstrumentindia.com

Description :

The Shock Absorber Test Stand trainer explores the construction and performance of shock absorbers as well as the operation and characteristics of suspensions. Shock Absorber Test Stand trainer is a modern up-to-date demonstration system and test stand for suspension and shock absorbers featuring a computer for data acquisition of test data and analysis of results.

Technical Specification :

Shock Absorber Test Stand



Laboratory Instrument India