



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Diesel Sleeve Metering Pump Cutaway

**Product Code :**  
LBNY-0005-120000218



**Description :**

The cutaway shows the internal components of an electronic control and cold advance. The Diesel Sleeve Metering Pump Cutaway provides a hands-on experience.

**Features:**

Hand Crank  
Powder Coated Metal Frame  
Late Model Recycled Sleeve Metering Pump - Caterpillar  
Diamond Plate Top Surface.

**Technical Specification :**

Diesel Sleeve Metering Pump Cutaway



**Laboratory Instrument India**