



Email : sales@laboratoryinstrumentindia.com

Product Name :
Diesel Sleeve Metering Pump Cutaway

Product Code :
LBNY-0005-120000218



Description :

The cutaway shows the internal components of an electronic control and cold advance. The Diesel Sleeve Metering Pump Cutaway provides a hands-on experience.

Features:

Hand Crank
Powder Coated Metal Frame
Late Model Recycled Sleeve Metering Pump - Caterpillar
Diamond Plate Top Surface.

Technical Specification :

Diesel Sleeve Metering Pump Cutaway



Laboratory Instrument India