



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hybrid Cutaway Trainer

Product Code :
LBNY-0005-120000220



www.laboratoryinstrumentindia.com

Description :

This training unit includes an electric motor which rotates the engine at a low RPM to show the operation and interaction of the internal working components. The Hybrid Cutaway Trainer is a complete Toyota Hybrid I-4 engine and planetary gear with sections cut away to expose the internal mechanical and electrical operating systems. A Reference Package is provided which includes service information specific to this engine. Intake, Exhaust, Cooling, and Lubrication systems are color coded to provide visual cues and allow for easy identification.

Technical Specification :

Hybrid Cutaway Trainer



Laboratory Instrument India