



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Hybrid Cutaway Trainer

**Product Code :**  
LBNY-0005-120000220



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

This training unit includes an electric motor which rotates the engine at a low RPM to show the operation and interaction of the internal working components. The Hybrid Cutaway Trainer is a complete Toyota Hybrid I-4 engine and planetary gear with sections cut away to expose the internal mechanical and electrical operating systems. A Reference Package is provided which includes service information specific to this engine. Intake, Exhaust, Cooling, and Lubrication systems are color coded to provide visual cues and allow for easy identification.

**Technical Specification :**

Hybrid Cutaway Trainer



**Laboratory Instrument India**