




Email : sales@laboratoryinstrumentindia.com

Product Name :
Electrical Machine Trainer Electronic Engineering Lab

Product Code :
LBNY-0002-620037



Description :

Electrical Machine Trainer Electronic Engineering Lab

Technical Specification :

Electrical Machine Trainer

Electric Machine, DC, AC, Shunt, Squarrel Cage,
Compound, Techo meter, Brake.

All machines are mounted on sturdy base frame with easy machine interchangeability

Use of gearless coupling facilitates screwless coupling.

With due emphasis on student safety, machines operate up to 300 W power levels and up to 1500 RPM without compromising on didactic use.

Easy and safe wiring by students using 4 mm sturdy shrouded banana patch cords and shrouded socket arrangements.

Trunnion mounted DC integrated machine is used as dynamometer for loading other machines (motors generators) unlike magnetic power brake or eddy current brake which can load only coupled motors and not generators.

Facility to measure shaft power using electronic torque / speed measurement.

Electrical Machine Trainer Electronic Engineering Lab, Electrical Machine Trainer Electronic Engineering Lab Equipments, Electrical Machine Trainer Electronic Engineering Lab Tools, Electrical Machine Trainer Electronic Engineering Lab Tool Kits, Electrical Machine Trainer Electronic Engineering Lab Manufacturers, Electrical Machine Trainer Electronic Engineering Lab Suppliers from India, China, Kenya



Laboratory Instrument India