




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Fuel Pump System

**Product Code :**  
LBNY-0005-120000231



**Description :**

**Features:**

- Easy Fuel Injection Volume Measurement
- A 15 Second Timer is Incorporated for Consistent Measurements at various RPM and Fuel Pressure Levels
- The Injector/Fuel Pump Trainer Demonstrates the Realistic Operation of Both a Fuel Pump Delivery System and a Fuel Injector
- DVOM/Oscilloscope Test Points for the Fuel Injector, Fuel Pump Terminals, and the Fuel Pump Relay Control Circuit
- Internal 12V Power Supply
- Adjustable Controls for Both RPM and Load Control
- Has an Adjustable Fuel Pressure Regulator (with Pressure Gauge) for
- Uses Actual Automotive Components Including the Fuel Pump, Fuel Injector, and Fuel Pump Relay
- Speed Variable from 200 - 5,140 RPM
- Courseware is Included (Student Manual and Instructor Guide).

**Technical Specification :**

Fuel Pump System

---

# Laboratory Instrument India

