



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Transfer Case Cutaway Trainer

**Product Code :**  
LBNY-0005-120000232



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

It provides hands-on experience to complete Tasks. The Transfer Case Cutaway Trainer presents the viewing and operation of a Transfer Case.

**Features:**

Handle allows rotation of gears  
Actual late-model transfer case  
Case cutaway painted to enhance viewing  
Mounted on a powder-coated metal stand.

**Technical Specification :**

Transfer Case Cutaway Trainer



**Laboratory Instrument India**