



Email : sales@laboratoryinstrumentindia.com

Product Name :
Transfer Case Cutaway Trainer

Product Code :
LBNY-0005-120000232



www.laboratoryinstrumentindia.com

Description :

It provides hands-on experience to complete Tasks. The Transfer Case Cutaway Trainer presents the viewing and operation of a Transfer Case.

Features:

Handle allows rotation of gears
Actual late-model transfer case
Case cutaway painted to enhance viewing
Mounted on a powder-coated metal stand.

Technical Specification :

Transfer Case Cutaway Trainer



Laboratory Instrument India