




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Fan Speed Control Trainer

**Product Code :**  
LBNY-0005-120000234



**Description :**

The Fan Speed Control Trainer circuits use the same wire colors as the manufacturers wiring diagram. Students follow printed schematics to wire up either a single-speed single fan system or a two-speed two fan system. Incandescent bulbs are substituted in place of the actual cooling fan motors for safety purposes. Fan Speed Control Trainer presents the hands-on construction, study and live operation of a modern Fan Speed Control System as found on late-model General Motors vehicles.

**Technical Specification :**

Fan Speed Control Trainer



**Laboratory Instrument India**