




Email : sales@laboratoryinstrumentindia.com

Product Name :
Power Seat Trainer

Product Code :
LBNY-0005-120000237



Description :

Power Seat Trainer allows for instructor demonstration, student hands-on reinforcement and troubleshooting practice. The Power Seat Trainer is part of an electrical systems program which presents the live operation and study of power seat systems. Powered by a 12VDC automotive battery or 12VDC/23A power supply.

Features:

- Intermittent and hard fault codes can be inserted using the built-in keypad, a single computer, or the Ability to perform actual service manual test procedures.
- Provides an advanced level of instruction by duplicating actual on-vehicle troubleshooting procedures.
- Courseware includes Instructor Guide, Student Manual, and Service Manual Information.
- Actual new vehicle components, connectors, and factory wire colors.
- Compatible with Instructor Management Program (IMP).
- Battery posts are protected against polarity reversal.
- Includes internal power supply for fault activation.

Technical Specification :

Power Seat Trainer

Laboratory Instrument India

