




Email : sales@laboratoryinstrumentindia.com

Product Name :
Diesel System Trainer

Product Code :
LBNY-0005-120000238



Description :

The trainer enables actual service test equipment to be used for diagnostic procedures, such as digital meters, breakout boxes, oscilloscopes and standard scan tools, including those having bi-directional and snap shot functions. The procedures used to diagnose trainer faults to be the same as performed on a vehicle. The program includes courseware for use by students and instructors. The Diesel System Trainer presents the real-world operation and study of the diesel fuel injection system.

Features:

- Electronic Control Module (ECM)
- Glow Plug Relay with Indicator Light
- Fuel Pump Relay with Indicator Light
- Engine Coolant Temperature Sensor
- Malfunction Indicator Light (MIL)
- Intake Air Temperature Sensor
- Fuse Block
- Standard Diagnostic Tool Connector (DLC)
- Service Throttle Soon Light.

Technical Specification :

Diesel System Trainer



Laboratory Instrument India