




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Diesel System Trainer

**Product Code :**  
LBNY-0005-120000238



#### **Description :**

The trainer enables actual service test equipment to be used for diagnostic procedures, such as digital meters, breakout boxes, oscilloscopes and standard scan tools, including those having bi-directional and snap shot functions. The procedures used to diagnose trainer faults to be the same as performed on a vehicle. The program includes courseware for use by students and instructors. The Diesel System Trainer presents the real-world operation and study of the diesel fuel injection system.

#### **Features:**

- Electronic Control Module (ECM)
- Glow Plug Relay with Indicator Light
- Fuel Pump Relay with Indicator Light
- Engine Coolant Temperature Sensor
- Malfunction Indicator Light (MIL)
- Intake Air Temperature Sensor
- Fuse Block
- Standard Diagnostic Tool Connector (DLC)
- Service Throttle Soon Light.

#### **Technical Specification :**

Diesel System Trainer



**Laboratory Instrument India**