




Email : sales@laboratoryinstrumentindia.com

Product Name :
Anti Lock Brakes Traction Control System

Product Code :
LBNY-0005-120000240



Description :

The ABS / TCS System trainer runs actual Service Manual test procedures, and all wire colors are consistent with the original vehicle schematics. This ABS / TCS System trainer provides for an advanced level of instruction by duplicating actual on vehicle operation and troubleshooting.

Features:

- Antilock Mode or Traction Control Mode Operation
- Class 2 "Bus" for Communications with Scan Tools
- Realistic Accelerator and Brake Pedal Operation
- Fault Insertion Locally with Keypad or Remotely
- Torque Request/Torque Delivered Indicators
- BPMV Solenoid Operation Indicators
- Selection of "Controlled Wheel" (LF, RF, LR, RR)
- Road Condition Selectable: Dry, Wet, Snow or Ice
- Courseware Included (Student Manual, Instructor Guide, Service Manual Information).

Technical Specification :

Anti Lock Brakes Traction Control System

Laboratory Instrument India

