



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Supplemental Inflatable Restraint Dual Air Bag System

**Product Code :**  
LBNY-0005-120000241



#### **Description :**

The Air Bag portion of the trainer is a nonfunctional disarmed unit. The Supplemental Inflatable Restraint Dual Air Bag System trainer uses an actual SIR System that electrically responds and operates as in the automobile.

#### **Features:**

- Provides an advanced level of instruction by duplicating actual on-vehicle troubleshooting procedures.
- Actual new vehicle components, connectors, and factory wire colors.
- Courseware includes Instructor Guide, Student Manual, and Service Manual Information.
- Crash simulation using sensor impact push switches and bag deployment strobe.
- Intermittent and hard fault codes can be inserted using the built-in keypad, a single computer.
- Ability to perform actual service manual test procedures.
- Performing airbag disposal procedures
- Study of removal and installation of system components.

#### **Technical Specification :**

Supplemental Inflatable Restraint Dual Air Bag System



**Laboratory Instrument India**

