




Email : sales@laboratoryinstrumentindia.com

Product Name :
Engine Repair Trainer

Product Code :
LBNY-0005-120000243



Description :

The Engine Repair Trainer provides an advanced level of instruction by duplicating actual on-vehicle service manual procedures. The Engine Repair trainer presents the rebuilding and service of automobile engines. The Engine Repair Trainer allows a 360 degree rotation of the engine.

Features:

- Mounted on a Rebuild Stand with Rolling Casters
- Late Model Recycled Engine Assembly
- Metal Pan Captures Oil Leaks
- Power Steering Pump
- Starter
- Flywheel
- Exhaust Manifold Assembly
- Safety Guard
- Alternator
- Injector Assembly
- All electronic components necessary for proper engine operation.

Technical Specification :

Engine Repair Trainer



Laboratory Instrument India