




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Engine Repair Trainer

**Product Code :**  
LBNY-0005-120000243



**Description :**

The Engine Repair Trainer provides an advanced level of instruction by duplicating actual on-vehicle service manual procedures. The Engine Repair trainer presents the rebuilding and service of automobile engines. The Engine Repair Trainer allows a 360 degree rotation of the engine.

**Features:**

- Mounted on a Rebuild Stand with Rolling Casters
- Late Model Recycled Engine Assembly
- Metal Pan Captures Oil Leaks
- Power Steering Pump
- Starter
- Flywheel
- Exhaust Manifold Assembly
- Safety Guard
- Alternator
- Injector Assembly
- All electronic components necessary for proper engine operation.

**Technical Specification :**

Engine Repair Trainer



**Laboratory Instrument India**