




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Supplemental Inflatable Restraint System

**Product Code :**  
LBNY-0005-120000247



### Description :

The Supplemental Inflatable Restraint System is based on a CAN platform vehicle. The Air Bag portion of the trainer is a non-functional disarmed unit. The Supplemental Inflatable Restraint System uses an actual SIR System that electrically responds and operates as in the automobile.

### Features:

- Intermittent and hard fault codes can be inserted using the built-in keypad, a single computer.
- Ability to perform actual service manual test procedures.
- New vehicle components, connectors, and factory wire colors.
- Provides an advanced level of instruction by duplicating actual on-vehicle troubleshooting procedures.
- Study of removal and installation of system components.
- Performing airbag disposal procedures.
- Includes internal power supply.
- Courseware includes Instructor Guide, Student Manual, and Service Manual Information.
- Keypad displays built-in scan tool data.
- Crash simulation using sensor impact push switches and bag deployment strobe.

### Technical Specification :

Supplemental Inflatable Restraint System

---

# Laboratory Instrument India

