




Email : sales@laboratoryinstrumentindia.com

Product Name :
Engine Cooling Fan System

Product Code :
LBNY-0005-120000248



Description :

Features:

- Capable of actual service manual test procedures
- Connects to scan tool from DLC connector
- Presents the live operation and study of a late-model System
- New-model vehicle components and connectors
- Faults can be inserted via keypad, local computer
- Courseware Included (Student Manual, Instructor Guide, Service Manual Information)
- Includes internal power supply for fault activation
- Power provided by 12 VDC battery or external power supply (option available).

Technical Specification :

Engine Cooling Fan System



Laboratory Instrument India