




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Engine Cooling Fan System

**Product Code :**  
LBNY-0005-120000248



#### **Description :**

#### **Features:**

- Capable of actual service manual test procedures
- Connects to scan tool from DLC connector
- Presents the live operation and study of a late-model System
- New-model vehicle components and connectors
- Faults can be inserted via keypad, local computer
- Courseware Included (Student Manual, Instructor Guide, Service Manual Information)
- Includes internal power supply for fault activation
- Power provided by 12 VDC battery or external power supply (option available).

#### **Technical Specification :**

Engine Cooling Fan System



**Laboratory Instrument India**