




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Diesel Injection Pump System Trainer

**Product Code :**  
LBNY-0005-120000249



#### **Description :**

The program provides an operational study of the diesel fuel injection system. The program includes courseware for use by students and instructors. This Diesel Injection Pump System trainer presents the real-world operation and study of the Stan dyne rotary diesel fuel injection pump.

#### **Features:**

- Crank Sensor Reluctor Wheel
- Crank Sensor with Terminal Test Points
- Diesel Fuel Injectors
- Lift Fuel Pump
- Stanadyne Rotary Injection Pump with:
  - a. Injection Control Module with Test Points
  - b. Optical Sensor with Terminal Test Points
  - c. Engine Shutoff Solenoid with Test Points
  - d. Injection Timing Motor with Test Points.

#### **Technical Specification :**

Diesel Injection Pump System Trainer

---

# Laboratory Instrument India

