



Email : sales@laboratoryinstrumentindia.com

Product Name :
Diesel Technology Trainer

Product Code :
LBNY-0005-120000250



Description :

It provides an advanced level of instruction by duplicating actual on-vehicle troubleshooting procedures. The Diesel Technology Trainer presents the live operation and study of a late model diesel pickup truck. Requires 12VDC automotive batteries.

Features:

- Diesel Particulate Filter System
- Portable via rolling casters
- Urea SCR Diesel Exhaust Fluid
- Compatible with actual service manual test procedures
- Performs DPF regeneration with scan tool
- Diesel Oxidation Catalyst
- Has a fully functional engine, transmission, front brake system, steering system, suspension system, fuel injection system, glow plug system, emission control, instrument panel, front and rear lighting system, and an air conditioning system
- Simulator/Fault Insertion Box for inserting realistic faults into the system

Technical Specification :

Diesel Technology Trainer

Laboratory Instrument India

