



Email : sales@laboratoryinstrumentindia.com

Product Name :
Automotive Technology Trainer

Product Code :
LBNY-0005-120000251



Description :

Features:

Assists in the Study and Diagnosis of most Automobile Systems.

The Automotive Technology Trainer or car clip is a late model recycled front-wheel-drive vehicle.

Including:

- a. Engine Performance
- b. Heating/Air Conditioning
- c. Suspension/Steering
- d. Electrical/Electronic
- e. Brake
- f. Transmission

Performs Actual Service Manual Diagnostic Procedures

Aids in Developing Real-World System Skills and Knowledge

Includes a Fault/Simulator Panel with Multiple Faulting Options and Features Including:

- a. Real-Time voltage readouts for Various Selectable Sensors
- b. Investigates Downstream Oxygen Sensors Control of Fuel Trim
- c. Voltage Adjustment Potentiometers for each Sensor
- d. Individual or Multiple Sensors are selectable for Voltage Substitution
- e. Demonstrates Catalytic Converter Punch-Through
- f. Interactive Demonstration/Study of Fuel Control by Oxygen Sensor Feedback.

Technical Specification :



Laboratory Instrument India