



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electronic Automatic Transmission Rebuild

Product Code :
LBNY-0005-120000257



Description :

Electronic Automatic Transmission Rebuild trainer presents the rebuilding and service of Electronic Transmissions.

Features:

Stand Includes Component Storage and Drain Pan
Powder-Coated Frame Mounted on Casters
Recycled Transmission
Rebuild Stand allows 360 Rotation of Transaxle
Service Manual Available.

Technical Specification :

Electronic Automatic Transmission Rebuild



Laboratory Instrument India