



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Electronic Automatic Transmission Rebuild

**Product Code :**  
LBNY-0005-120000257



**Description :**

Electronic Automatic Transmission Rebuild trainer presents the rebuilding and service of Electronic Transmissions.

**Features:**

Stand Includes Component Storage and Drain Pan  
Powder-Coated Frame Mounted on Casters  
Recycled Transmission  
Rebuild Stand allows 360 Rotation of Transaxle  
Service Manual Available.

**Technical Specification :**

Electronic Automatic Transmission Rebuild



**Laboratory Instrument India**