



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electricity and Electronics Fundamentals Trainer Electronic Engineering Lab

Product Code :
LBNY-0002-620040



Description :

Electricity and Electronics Fundamentals Trainer Electronic Engineering Lab

Technical Specification :

Electricity and Electronics Fundamentals Trainer

Electricity and Electronics Fundamentals Trainer covering the concepts, theory and applications of basic electricity and electronics that are required for technicians. This comprehensive course is presented in a systematic way starting the student out at the beginning covering basic theory and measurements with meters.

Trainer board, component and Circuit based PCB Module

A. Base Station: Fixed regulated DC Power Supplies: +5V, +15V, -15V

Variable Regulated DC Power Supply 0 - 25V

Potentiometers: 4 units

Pulse Switches: 2 units

8 bits LED display.

Toggle Switches: 5 spdt

Speaker: 2.25 inch, 8 ohm/0.25w

One set of 2 mm stackable test leads to perform all required inter-connections.

B. Set of PCB Board Experiment Modules.

Electricity and Electronics Fundamentals Trainer Electronic Engineering Lab, Electricity and Electronics Fundamentals Trainer Electronic Engineering Lab Equipments, Electricity and Electronics Fundamentals Trainer Electronic Engineering Lab Tools, Electricity and Electronics Fundamentals Trainer Electronic Engineering Lab

Tool Kits, Electricity and Electronics Fundamentals Trainer Electronic Engineering Lab Manufacturers, Electricity and Electronics Fundamentals Trainer Electronic Engineering Lab Suppliers from India, China, Kenya



Laboratory Instrument India