



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electronic Automatic Transaxle Rebuild

Product Code :
LBNY-0005-120000260



Description :

Electronic Automatic Transaxle Rebuild trainer presents the rebuilding and service of Electronic Transaxles.

Features:

- Rebuild Stand allows 360 Rotation of Transaxle
- Powder-Coated Frame Mounted on Casters
- Recycled Transaxle
- Stand Includes Component Storage and Drain Pan
- Service Manual Available.

Technical Specification :

Electronic Automatic Transaxle Rebuild



Laboratory Instrument India