




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Power Seat System

**Product Code :**  
LBNY-0005-120000268



### Description :

Power Seat System trainer allows for instructor demonstration, student hands-on reinforcement and troubleshooting practice. The Power Seat System is part of an electrical system program which presents the live operation and study of power seat systems. The Power Seat System is based on a CAN platform vehicle. Powered by a 12VDC automotive battery or 12VDC/30A power supply (not included).

### Features:

- Provides an advanced level of instruction by duplicating actual on-vehicle troubleshooting procedures.
- Ability to perform actual service manual test procedures.
- Compatible with Instructor Management Program (IMP).
- Courseware includes Instructor Guide, Student Manual, and Service Manual Information.
- Intermittent and hard fault codes can be inserted using the Instructor Management Program (IMP) via a computer connected to the trainer.
- New vehicle components, connectors, and factory wire colors.
- Hard faults can be inserted by using the built in keypad on the Power Seat System.
- Includes external power supply for fault activation.
- Battery posts are protected against polarity reversal.

### Technical Specification :

Power Seat System



**Laboratory Instrument India**