



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Digital Logic Trainer Electronic Engineering Lab

**Product Code :**  
LBNY-0002-620041



**Description :**

Digital Logic Trainer Electronic Engineering Lab

**Technical Specification :**

Digital Logic Trainer

Digital Logic Trainer (TTL) Logic Trainer Board For Electrical Lab Training Digital Logic Trainer (TTL) /Logic Trainer Board based on 74 series has been designed specifically to make the students familiar with the study of TTL ICs and verification of the truth table of logic gates, flip-flops, Gated & Master Slave JK flip-flops, Schmitt Trigger, Expanders, Binary address, Counters, Shift registers, Multiplexer (Encoder), Demultiplexer (Decoder), 8Bit D/A Converter and 8Bit A/D Converter etc.

Basic and Advance Modules

Fixed regulated DC power supplies : +5V, +15V, -15V, variable 1.5 ~ 15VDC

Logic indicators : 16 logic indicator

Data switches : 16 (for TTL & CMOS)

Toggle switches : 4 (SPDT)

Push buttons : 3 (2NO, 1NC)

Pulsers : 2 pulsers (A,A and B,B)

Seven segments : 4 seven segment digital displays

Clock generator : 1Hz 100kHz

Test leads : 1 set of stackable 2mm test leads

Digital Logic Trainer Electronic Engineering Lab, Digital Logic Trainer Electronic Engineering Lab Equipments, Digital Logic Trainer Electronic Engineering Lab Tools, Digital Logic Trainer Electronic Engineering Lab Tool

---

Kits, Digital Logic Trainer Electronic Engineering Lab Manufacturers, Digital Logic Trainer Electronic Engineering Lab Suppliers from India, China, Kenya



## Laboratory Instrument India