



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
ATF Changer

**Product Code :**  
LBNY-0005-1210005



**Description :**

The ATF Changer equipment can be convenient, fast and efficient completing circulation Cleaning and fluid equivalent exchanging for old and new oil.

**Features:**

The usage of integrated functional pipeline modules ensures reliable operation of equipment  
A variety of alarm/ indication information and automatic protection measures prevents damage to Automatic Transmission System.

The electronic valve controls the flow direction, preventing cross contamination of the new and old fluid.

Monitored by electronic load cell and flow sensor, to double ensure that replacement precision.

Artistic Shape, easy to move.

Power supply has positive and negative self correcting recognition function, to avoid damage to equipment control unit.

**Technical Specification :**

ATF Changer



**Laboratory Instrument India**

