



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wheel Aligner

Product Code :
LBNY-0005-1210008



Description :

The Wheel Aligner equipment is used for detecting the alignment parameters of the automobile wheel, and comparing with the original design parameters, guiding the users to adjust the wheel alignment parameters accordingly.

Features:

Adopt Inclination sensor, Bluetooth wireless communication
High accuracy and high speed
Easy operation and learn.

Technical Specification :

Measurement Option/Measuring Accuracy

Front Camber-5

Rear Camber-5

Front toe-10

Rear toe-10

Include Angle-20

Setback-5

Castor-20

Thrust-1.5

Temperature:(10-45)

Total measuring range:(2'-4')

Dimensions

Length-800mm (10%)
Width-670mm (10%)
Height-1290mm (10%).



Laboratory Instrument India