



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Double Column Lifting Machine

**Product Code :**  
LBNY-0005-1210009



**Description :**

Wire rope equalization system forces synchronous movement of two carriages and effectively prevents vehicle from tilting.

Double hydraulic cylinder and high strength chain drive makes stable and smooth lifting and lowering.

Door rubber pad protects car door from being damaged in collision.

More lifting pads with different height meets the requirements for light lorries and SUV.

24V low voltage control box ensures operator's safety.

24V cylinder travel limit switch effectively protects cylinder and extends its lifespan.

**Technical Specification :**

Double Column Lifting Machine



**Laboratory Instrument India**