



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tyre Changer

Product Code :
LBNY-0005-12100011



Description :

It makes the automobile maintenance process more convenient and smooth. It's an automobile maintenance equipment using for auxiliary disassembly and assembly for automobile tires in vehicle maintenance.

Features:

Vertical and horizontal arms with pneumatic dual locking device.
Tool box under the pole
Slide with step function
Self-centering device.
Pole is tilt-able.

Technical Specification :

Rim clamping from outside: 12"-22"
Rim clamping from inside: 14"-24"
Motor Power: 1.1 /0.75 kw
Power Supply: 110/220/380V-3ph
Max Wheel diameter: 960 mm
Operating Pressure: 8-12 bar
Bead-breaker force: 2500kg

Laboratory Instrument India

