



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Tyre Changer

**Product Code :**  
LBNY-0005-12100011



**Description :**

It makes the automobile maintenance process more convenient and smooth. It's an automobile maintenance equipment using for auxiliary disassembly and assembly for automobile tires in vehicle maintenance.

**Features:**

Vertical and horizontal arms with pneumatic dual locking device.  
Tool box under the pole  
Slide with step function  
Self-centering device.  
Pole is tilt-able.

**Technical Specification :**

Rim clamping from outside: 12"-22"  
Rim clamping from inside: 14"-24"  
Motor Power: 1.1 /0.75 kw  
Power Supply: 110/220/380V-3ph  
Max Wheel diameter: 960 mm  
Operating Pressure: 8-12 bar  
Bead-breaker force: 2500kg

---

# Laboratory Instrument India

