



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Digital Circuit Training Kit

**Product Code :**  
LBNY-0005-12100015



**Description :**

Digital Circuit Training Kit it uses high level quality components, it is designed for carrying out multi-type digital circuit experiments and it is suitable for electronic technology research and development in the vocational and technical school and universities.

**Technical Specification :**

Digital Circuit Training Kit



**Laboratory Instrument India**