



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mobile Communication Training System

Product Code :
LBNY-0005-12100016



Description :

Output power : +5 to +33 db, 3.2 min to 2 w

Frequency bank: 935 to 960 MHZ

Cellular system: GSM receive

Power(trainer): 230V AC

Dimension (trainer) :41.5* 31.5 *13.5

Weight(trainer): 5kg

Accessories: loud speaker box, connection strip (mobile connector cable), mobile hand set, charger cord.

Mains housing (trainer): It is housed in an elegant cabinet.

Technical Specification :

Mobile Communication Training System



Laboratory Instrument India