



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mobile Communication Training System

**Product Code :**  
LBNY-0005-12100016



**Description :**

Output power : +5 to +33 db, 3.2 min to 2 w

Frequency bank: 935 to 960 MHZ

Cellular system: GSM receive

Power(trainer): 230V AC

Dimension (trainer) :41.5\* 31.5 \*13.5

Weight(trainer): 5kg

Accessories: loud speaker box, connection strip (mobile connector cable), mobile hand set, charger cord.

Mains housing (trainer): It is housed in an elegant cabinet.

**Technical Specification :**

Mobile Communication Training System



**Laboratory Instrument India**