



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cabinet Air Conditioner Training Device

**Product Code :**  
LBNY-0005-12100018



**Description :**

The Cabinet Air Conditioner Training Device can be used to train students to master the theoretical knowledge and technical skills of air conditioning installation, commissioning and maintenance.

**Technical Specification :**

Rated power: 1.5KW  
Input power: AC220V, 50HZ  
Dimension: about 950\*800\*1750mm



**Laboratory Instrument India**