



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Portable Recovery Unit

**Product Code :**  
LBNY-0005-12100019



**Description :**

**Features:**

Suitable for HCFC and HFC refrigerants and blends  
Spare parts for 1 year including 5 replaceable filters  
Overpressure protection system  
3/4 HP compressor  
Compact and easy to carry  
OFP feature  
Replaceable filter  
Size: 465\*225\*360mm

**Technical Specification :**

Portable Recovery Unit



**Laboratory Instrument India**