



Email : sales@laboratoryinstrumentindia.com

Product Name :

Resistance Decade Box Electronic Engineering Lab

Product Code :

LBNY-0002-620043



Description :

Resistance Decade Box Electronic Engineering Lab

Technical Specification :

Resistance Decade Box

Variable from 1 to 10 m Ω Resistance

The Decade Resistance boxes allow quick and easy setting of a wide range of resistance.

The main function of these decade boxes is to create resistance with a specific value by using a combination of rotary decade switches.

Resistance boxes are used in laboratory testing and calibration.

Resistance Decade Box Electronic Engineering Lab, Resistance Decade Box Electronic Engineering Lab Equipments, Resistance Decade Box Electronic Engineering Lab Tools, Resistance Decade Box Electronic Engineering Lab Tool Kits, Resistance Decade Box Electronic Engineering Lab Manufacturers, Resistance Decade Box Electronic Engineering Lab Suppliers from India, China, Kenya



Laboratory Instrument India

