



Email : sales@laboratoryinstrumentindia.com

Product Name :
Nitrogen Regulator

Product Code :
LBNY-0005-12100021



Description :

Features:

Material: Copper
Measurement range: (0-25) MPa
Inlet pressure: 15MPa
Working pressure 40 bar
Input thread: mm "G5/8
Output thread: M16X 1.5mm
Outlet pressure: 4MPa
Regulating range: 0.4-4MPa
Rated flow: 16m³/h
Weight: 2kg.

Technical Specification :

Nitrogen Regulator



Laboratory Instrument India

